eBike L & EL

Semi-Couch & Stress Echo Ergometer for load echo cardiography

The new eBike L & eBike EL are specially designed for couch ergometry and stress echo applications. Both ergometers have a wide range of uses. The eBike L & eBike EL have a very ergonomic design that assures maximum patient safety and comfort. The proven eBike measurement intelligence combined with a very easy to use operator interface convinces on the clinical side and during the daily use. Both eBike L & eBike EL are offered in different combinations. It is offered with or without integrated Blood Pressure module and with standard or separate Control Terminal. The integrated Blood Pressure module can also be offered as an upgrade.

The eBike L & eBike EL has an automatic electrical adjustment device, which enables the ergometer to be quickly moved from the flat position directly to a position of 45° (also lateral 0° to 45° for eBike EL). When the patient is in this semirecumbent position, the muscular regions and the arms are relaxed, thus providing the ideal conditions for taking ECG and blood pressure measurements. Should the need arise during an examination, the ergometer, with its patient, can be converted from the semi-recumbent to the flat position within seconds, enabling treatment to be carried out immediately.





The hinged couch is mounted on a firm base. There are two footplates to facilitate access for elderly or disabled patients. The legs are positioned at the ideal angle to the tread crank and there is no danger of the feet slipping. The headrest, complete with shoulder supports and arm rests, is secured by a guide rail and this ensures that blood pressure measurement can be taken safely. The computer-controlled blood pressure module ensures maximum accuracy and includes automatic failsafe measurements.

The examination of a patient under stress is increasingly asked for in the special line of echo cardiography as well as in electrocardiography in order to render a reliable diagnosis possible. The leaning/halfleaning safety ergometer eBike EL from GE Healthcare has been designed especially for echo cardiography.



Specifications eBike L & eBike EL

Operating Mode Continuous operation
Power Supply 110 to 120 V, 50 to 60 Hz

220 to 240 V, 50 to 60 Hz

Power Consumption 80 VA max.

Braking Principle Computer-Controlled eddy current

bracke with torque measurement;

speed independent to DIN VDE 0750-238

Load Range 20 to 999 Watt,

speed (RPM)-independent

RPM Range 30 to 130 RPM

Load Error ±5% between 20 and 900 W and

±3W max. between 20 and 100 W

Load Increments Configurable: 1, 5, 10 or 25 W

Moment of Inertia10 kg x m²Flywheel Mass7 kgCrank Length170 mm

Internal Protocols 5 fixed protocols, 10 user-configu-

rable protocols

Interfaces 2 x RS232: 9-pin Sub-D

1 x RS232: 5-pin DIN socket Analog input for target load:

8-pin DIN socket

Analog output for current load:

8-pin DIN socket Remote start EKG unit: 1...30 s before load change 115 x 88 mm, 320 x 240 pixels

Environment/Operation

Display

Temperature: between +10°C and +40°C Relative Humidity: 30 to 75%, no condensation

Atmospheric Pressure: 700 to 1060 hPa

Certification CE 0123

Accessories

Please contact your GE representative for our wide range of tested and approved quality accessories and

consumables.

eBike L

Patient Weight 160 kg max.

Saddle Height Adjustment

Motor-driven, continuous for patients

from 120 to 210 cm

Tilt Range Motor-driven, from a flat position to 45°

Dimensions (WxL): 800 x 2350 mm min. (tilted 45°,

head support in bottom position) 800 x 2520 mm max. (horizontal, head support in upper position)

Weight approx. 110 kg

eBike EL

Patient Weight 140 kg max.

Saddle Height Adjustment

Motor-driven, continuous for patients

from 120 to 200 cm

Tilt Range Motor-driven, from a flat position to

45°; lateral 0° to 45°

Dimensions (WxL): 800 x 2100 mm min. (tilted 45°,

head support in bottom position) 1200 x 2600 mm max. (tilted 45°, head support in upper position

Weight approx. 140 kg

GE Medical Systems Information Technologies GmbH Munzinger Strasse 3 79111 Freiburg, Germany Tel. +49 761 4543 0 • Fax +49 761 4543 233

www.gehealthcare.com

© 2005 General Electric Company – All rights reserved. GE and GE Monogram are trademarks of General Electric Company.

serieral Electric company.

GE Medical Systems Information Technologies GmbH, a General Electric company, going to market as



